## Instructions For Use

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Catalog Number	<b>803-S 407</b>	Effective Date	<b>August 12, 2021</b>	Status	<b>APPROVED</b>
Doc ID	<b>IFU-IMP-0006</b>	Version	<b>3</b>	Page 1 of 8	

**Intended Use/Intended Users:** The De Mayo Knee Positioner’s intended use is to position the patient’s leg during a surgical procedure, including but not limited to total knee replacement, partial knee replacement and revision. It is also intended to be used in Stryker Robotic Systems. Orthopedic surgeons are the intended users.

**Target Patient Group:** Patient selection depends on the judgment of the surgeon. Surgeon must consider size of the patient. The heel to popliteal region should be in between 13 inches (33 cm) and 18 inches (46 cm) and tibia should be over 11 inches (28 cm) in length.

**Contraindications:** This device is not designed, sold or intended for use except as indicated.

### Warnings/Precautions

- Do not use device in a manner that doesn’t follow these instructions for use.
- Untrained personnel - review and understand the IFU
- Incorrect set up – follow IFU
- Safety: Always use IMP Patient Protective Pads
- **DO NOT PLACE POSITIONER IN A MILK BATH OR LUBRICATE**
- Max number of reuses:
  - Unplanned carriage movement
  - Excessive medial lateral movement of boot when the carriage is locked
  - Unrepairable wear and/or damage to the baseplate including locking pins

### Risks:


- Device should not come in direct contact with the patient. Patient is protected with Patient Protective Pad
- Sterilized by end user
- No shelf life issues
- Re-usable
- Must be in compliance with the IFU with IMP accessories and components

### Complaints and Adverse Events:

For complaints and adverse events, contact IMP, EU representative, and the appropriate regulatory authorities for specific country.

Innovative Medical Products  
 87 Spring Lane  
 Plainville, CT 06062  
 P: 860-793-0391  
 F: 866-459-1805  
 info@impmedical.com

Advena Ltd  
 Tower Business Centre, 2nd Flr,  
 Tower Street, Swatar, BKR 4013 Malta

### Product Identification

Part No	Product Name	GTIN
803-S 407	De Mayo Knee Positioner <sup>®</sup> 25" w/ Composite Boot, Single Lever Clamp and Stationary Vertical Pins	00696588000787
M-408-10	MAKO Sterile Patient Protective Pads for Knee Positioners 10/case	00696588000916

**Consumables:**

- Sterile Patient Protective Pads
- Coflex

**Disposal of unit:**

If a device is being returned for evaluation, please contact your local MAKO representative for shipping/handling information. If the device is not being returned to MAKO Surgical Corp., instruments are to be disposed of in accordance with applicable laws, rules, and regulations for the disposal of biohazardous waste. Follow all guidelines for biohazardous waste in accordance with the Centers for Disease Control and Prevention guidelines as well as applicable federal/national, state and local regulations.

For additional information on repair or disposal, contact IMP at [sales@impmedical.com](mailto:sales@impmedical.com).

**Acceptable Accessories:**

Part No	Product Name	GTIN
907	De Mayo Universal Distractor®	00696588001388
713-717	De Mayo Push Button Clamp -Locking Pin	00696588001517
120	Quad™ Clamp	00696588006178
803-ABD	Aluminum Boot with Distractor Block	00696588001463
803-ABDA	Aluminum Boot with Distractor Block 110°	00696588006468
919	Friction Fold Down Pin Clamp Base	00696588005829

**Symbol Glossary**

Symbol	Title	Description	Standard
	Authorized Representative in the European Community	Indicates the authorized representative in the European Community.	ISO 15223-1:2016
	Batch Number	Indicates the manufacturer's batch code so that the batch or lot can be identified.	ISO 15223-1:2016
	Catalogue number	Indicates the manufacturer's catalogue number so that the medical device can be identified. The manufacturer's catalogue number shall be placed after or below the symbol and adjacent to it.	ISO 15223-1:2016
	Caution	Indicates the need for the user to consult the instructions for use for important cautionary information such as warnings and precautions that cannot be presented on the medical device itself.	ISO 15223-1:2016
	Complies with European Directives		
	Consult instructions for use		ISO 15223-1:2016
	Date of manufacture	The date must be presented in the following format: YYYY-MM-DD	FDA 21 CFR 801

Symbol	Title	Description	Standard
	Manufacturer	This symbol shall be accompanied by the name and address of the manufacturer.	ISO 15223-1:2016
	Medical device		ISO 15223-1:2016
	Non-sterile	This symbol should only be used to distinguish between identical or similar medical devices sold in both sterile and non-sterile conditions.	ISO 15223-1:2016
	Not made with natural rubber latex		Manufacturer defined
	Serial Number	The manufacturer's serial number shall be placed after or below the symbol and adjacent to it.	ISO 15223-1:2016
	Unique device identifier		ISO 15223-1:2016

1. Replace OR Table Pads with green <i>imp</i> <sup>®</sup> Pads. <b>Note: It is not necessary to remove X-Ray cassette plates.</b>	
2. Position the patient with the gluteal fold at the separation of the pads.	
3. During prep of the surgical leg, remove single OR Table Pad.	
4. Create a 3"-5" fold in the final drape under the buttocks.	
5. Place De Mayo Knee Positioner <sup>®</sup> in the well to achieve maximum flexion.	
6. Check flexion for final approval of freedom of the drape to allow full range of motion, when the positioner is locked on the drapes.	

7. Position Single Lever Clamp directly under guide pins and lock to the OR Table Side Rail utilizing either

**See Clamp IFU on page 4**

8. Insert imp® Patient Protective Pad into the sterile boot.

**For patient's safety, always use imp® Patient Protective Pads®.**

9. Wrap cohesive bandage around the foot starting with a minimum of six (6) foot wraps, tear, and finish the wrap above the distractor block around the calf.

**Note: Do not wrap cohesive bandage over the distractor block on the back of the boot.**

**After surgery - recommend spraying the positioner with a pre-treatment prior to resterilization**

## Safety Test

Tighten clockwise to lock the knob and lever

1. Check that all knobs and handles are locked
2. Patient is fully protected by imp® Patient Protective Pads and cohesive wrap
3. Patient's leg can reach full flexion
4. Check that both Guide Pins are fully seated


## IFU: De Mayo Single Lever Clamp


1. Position the clamp onto the O.R. table rail and with the lever in "unlock" position.

2. Rotate lever to "lock" position.

3. Squeeze jaws tight on side rail and tighten knurled knob clockwise all the way.

4. While holding clamp, rotate lever to “unlock”.


5. Turn knurled knob clockwise 1/4 to 1/2 rotation.

6. Rotate lever to “lock” position requiring only 2 finger pressure to achieve locking. Works like a vice grip and will “snap” into place.


7. Slide and guide pins into holes of the clamp.

## Safety Test

1. Ensure clamp knob is tightened
2. Lever in the “locked” position to the right side of the clamp
3. Check that both guide pins are fully seated



## Spring-Loaded Handle Re-Tightening Instructions:

If the boot and carriage require additional tightening after locking into place, follow the instructions below to re-tighten. **Do Not Remove Center Screw.**

1. Starting with the handle in the locked position as shown, but ball of boot and carriage are not tight enough,

2. <b>Pull</b> spring-loaded handle toward you		
3. <b>Rotate</b> , counter-clockwise, while holding in the pulled-out position		
4. <b>Release</b> , at desired position snapping into gear, and <b>Tighten</b> , until locked		

## Cleaning and Sterilization Procedure (hardware)

NOTE: ALL SOLUTIONS MUST BE COMPATIBLE WITH ALUMINUM & STAINLESS STEEL

### Recommended Washer / Decontaminator Instructions:

- Soak the product in an enzyme solution. Follow the manufacturer's direction for dilution and soaking time.
- Put through washer / decontaminator according to manufacturer's instructions **with a detergent up to a pH of 9.0**  
**NOTE:** Select cycle that does not include lubrication.

### Recommended Hand Cleaning Instructions:

- Pre-Soak the product / components in an enzyme solution. Follow the manufacturer's direction for dilution ratio and soaking time.
- Rinse the product in warm tap water.
- Wash the product with an instrument detergent up to a pH of **9.0** or enzyme product
- Rinse the product in warm tap water
- Soak or wipe the product down with a hospital approved and / or EPA approved germicide according to instructions.
- Rinse the product in warm tap water
- Dry thoroughly and wrap

### Recommended Sterilization Instructions:

- Ensure that all parts are thoroughly cleaned.
- Make sure that all movable parts are loose and can move freely.
- **DO NOT PLACE POSITIONER IN A MILK BATH OR LUBRICATE**
- Double wrap in two disposable wraps. Use 48 inch x 48 inch wraps. (Approximately 122 cm x 122 cm).
- Run normal vacuum cycle for your institution
- **STEAM STERILIZATION ONLY- ALL OTHER STERILIZATION METHODS NOT VALIDATED**

**MINIMUM PARAMETERS PRE-VAC STERILIZATION**

Product	Part Number	Temperature Setting	Exposure Time	Dry Time
De Mayo Knee Positioner	803-S-407	270°F to 272°F 132°C to 134° C	4 minutes	20 minutes
Single Lever Clamp	903	270°F to 272°F 132°C to 134° C	4 minutes	20 minutes

**CAUTION: The positioner must be cool before applying to the patient**

Sterilization Parameters Certified by:

- Micro Test Laboratories (now Accuratus Lab Services)
- Accuratus Lab Services
- HIGHPOWER Validation Testing and Lab Services

**Cleaning of the OR Table Pad**

Approved products and methods for cleaning and disinfecting of OR Table pads CAT# **398-LG** OR Table Pads and **398-ABG** Arm Board Pads sold by Innovative Medical Products.

**DO NOT APPLY TAPE TO ANY PAD****CLEANING OF SOILS & STAINS**

Wash with a generous application of neutral soap suds and lukewarm water up to 103 °F (39 °C). Rinse with water and dry. Do not immerse.

or

Use a diluted bleach solution of 10-parts water to 1-part bleach. Wipe dry with a clean cloth. Do not immerse. Bleach **MUST** be diluted to the 10:1 ratio.

All chemicals/solutions **MUST** be diluted to the manufacturer's recommendations. No chemical or solution should be left to dry on the pad surface. After the allotted time has passed to ensure disinfection, the pad surface should be wiped again with a wet cloth to remove the remaining chemical or solution.

**DISINFECTING**

Use Lysol Brand III I.C. Disinfectant Aerosol Spray by Reckitt Benckiser Inc., per manufacturer's instructions.

or

Use a diluted bleach solution of 10-parts water to 1-part bleach. Wipe dry with a clean cloth. Do not immerse. Bleach **MUST** be diluted to the 10:1 ratio.

All chemicals/solutions **MUST** be diluted to the manufacturer's recommendations. No chemical or solution should be left to dry on the pad surface. After the allotted time has passed to ensure disinfection, the pad surface should be wiped again with a wet cloth to remove the remaining chemical or solution.

**The following Steris products are also approved for use in cleaning and disinfecting our pads:**

- Steris Coverage Plus wipes 1608-WC, -GS, -G4
- Steris Coverage Spray HB Plus 1624-77
- Steris Coverage Spray TB Plus 1629-B4, -08

**The use of ANY cleaning/disinfecting chemical other than those on this list will void the warranty on any Table Pads sold by IMP.**

Scan for language specific IFUs

(EN, DA, KO, DE, EL, ES, FL, FR, ZH-CN, IT, NL, NO, PL, PT, PT-BR, RO, SV, TR)

(01) 0 0696588 00078 7  
(11) 210203  
(21) AAA

De Mayo Knee Positioner® is a registered trademark of Innovative Medical Products, Inc.  
imp® products are protected by patent & patent pending rights ~ go to [impmedical.com/patents](http://impmedical.com/patents)  
© August 2021 Innovative Medical Products, Inc.  
ALL RIGHTS RESERVED

**Innovative Medical Products**  
87 Spring Lane  
Plainville, CT 06062  
P: 860-793-0391  
F: 866-459-1805

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
Advena Ltd  
 Tower Business Centre, 2nd Flr,  
 Tower Street, Swatar, BKR 4013 Malte

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