

# imp De Mayo & V<sup>2</sup>E<sup>tm</sup> Knee Positioner Configurations

Versatility with V<sup>2</sup>E<sup>tm</sup>

<u>Description</u>	<u>Components</u>
V <sup>2</sup> E <sup>tm</sup>	Base Plate, Boot, Clamp (#114-CL), Extension
V <sup>2</sup> tm	Base Plate, Boot, Clamp (#114-CL)
D <sup>2</sup>	Base Plate (with removable carriage), Boot, Clamp (#903, 713, 302)
D	Base Plate, Boot, Clamp (#903, 713, 302)



V<sup>2</sup>E<sup>tm</sup>

V<sup>2</sup>tm

D<sup>2</sup>

D

### Specific Features

Extension (#114-EXT)	✓			
Removable Carriage	✓	✓	✓	
Lightweight Base Plate	✓	✓		
Varus tilt control lock	✓	✓		
180° Flip Handle	✓	✓		
Single Bar Clamp Lock	✓	✓		
Clamp (#114-CL)	✓	✓		
Clamp (#903, 302, 713)			✓	✓
Fold Down Pins			✓	✓
25" Base Plate	✓	✓		✓
30" Base Plate	✓	✓	✓	✓
Composite Boot, 110°	✓	✓	✓	✓
Aluminum Boot, 90°	✓	✓	✓	✓
Spring Handle	✓	✓	✓	✓
E-Brake	✓	✓	✓	✓

Note: Top rail of the "D<sup>2</sup>" is wider than the "V Versions". The carriages are NOT interchangeable.  
 The Removable Carriage of the "D<sup>2</sup>" does NOT fit on the "D" Base Plate.