



Innovative Medical Products’ family of hip positioners provides simple, adjustable positioning for any sized patient, on any OR table

IMP’s system saves time and lowers risk for every hip procedure, putting patients in a better position for positive outcomes

PLAINVILLE, Conn. (July 20, 2018) – Innovative Medical Products’ family of hip positioners gives surgeons and OR staff everything they need to position patients faster, more securely and with less strain on personnel.

Developed by Edward De Mayo, M.D., an orthopedic surgeon, the De Mayo Hip Positioner® system is based on a secure patient positioning device that enables OR personnel to “build” the positioner around the patient, reducing stress on OR staff. The IMP system saves time and lowers risk for every hip procedure.

The system’s collection of stabilizing bars and support pads offers the right option for any patient body type – even the most challenging of obese patients. With the IMP positioner, surgeons get unobstructed access to the surgical site with support and security for obese patients that can’t be found in any other device.

“From establishing the position and securing the patient, to performing procedures and testing range of motion, we created a system that’s faster and easier while also giving doctors completely unrestricted access to the surgical field,” notes Dr. Edward De Mayo.

Besides the De Mayo Hip Positioner®, IMP’s family of hip positioners includes:

The **Universal Lateral Positioner®** employs a versatile system of plates, arms, swivels and pads to maximize positioning options. Successfully positioning patients since 1984, IMP’s Gold Standard for lateral positioning has a complete choice of options to meet any patient need; plus, it fits into a convenient hard case for easy storage and transportation.

The **MorphBoard®** pegboard system is the most stable patient positioning pegboard available for lateral decubitus positions. Modular components allow simple conversion for obese patient support on any OR table and multiple peg lengths provide a versatile, secure fit. The IMP solution comes with Patient Protective Pads to protect from cross-contamination.

The **De Mayo Expand-A-Peg® Locking Peg** delivers a new level of pegboard-positioning safety. The locking mechanism prevents loosening and motion to keep a patient securely in place during surgery. Pegs lock with a simple twist. The pegs are compatible with IMP’s MorphBoard® and other surgical pegboards.

IMP's **SlimLine™ Abduction Pillows** improve patient comfort and simplify post-hip-surgery care. The proprietary convoluted foam structure is more comfortable and flexible for patients. Large straps secure easily and hold legs safely in place. One size fits most patients, lowering inventory costs.

Customers: for more information about IMP's Family of Hip Positioners, please contact Innovative Medical Products, 87 Spring Lane, P.O. Box 8028, Plainville, CT 06062; phone: 860-793-0391 or toll free: 800-467-4944; or visit www.impmedical.com.

About Innovative Medical Products

Innovative Medical Products Inc.® designs, manufactures, and brings to market patient positioning devices for healthcare and sports related surgeries. With more than 40 years experience in supporting hospitals, surgery centers, and clinics, IMP provides positioning products that stabilize and protect patients during navigated, computer assisted and robotic surgeries. IMP devices save hospital and surgery center costs by minimizing staff, decreasing liability and improving the overall patient experience. IMP manufactures, patents and trademarks products to Federal and International standards for quality and safety for both the patient and the healthcare institution. For more information, visit: www.impmedical.com

###

Media contact: Harry McBrien, Maier True Communication; 860-470-3127; harry@maier.com