



## **Innovative Medical Products' De Mayo Knee Positioner® is the 'Gold Standard' in Family of IMP Knee Positioners**

### *Positioner features 'Push Button Release Mechanism'*

PLAINVILLE, Conn., (January 23, 2018) – About 15 years ago, IMP began manufacturing what was to become the company's "gold standard" in patient positioning for knee surgeries – the De Mayo Knee Positioner®, a solution that allows for precise control of flexion, internal and external rotation, and extension of the knee during surgery. Since its introduction, the IMP positioner has generated unprecedented sales and sparked the beginning of a steady evolution in IMP's providing the latest in patient positioning technology and devices for knee surgery.

Used in conjunction with the De Mayo Universal Distractor®, the De Mayo Knee Positioner® is like having "another pair of hands" in the operating room, freeing up OR staff for other time-saving responsibilities.

The De Mayo Knee Positioner® is secured to the OR table by the De Mayo Push Button™ Clamp, an operating table rail clamp that significantly decreases OR set-up time. The clamp is fastened to the operating table simply by squeezing its top and bottom jaws closed on the table's rail and then locking into place with a twist of the device's knob. The clamp is quickly and easily removed by pushing the release buttons on the device, lowering the bottom jaw of the clamp.

Manufactured in two different configurations, the De Mayo Push Button™ Clamp-Standard can be used with the De Mayo Knee Positioner®, De Mayo Ankle Distractor™, or any of the IMP table accessories that require a dual post locking clamp for sterile applications. The De Mayo Pin Locking Push Button™ Clamp adds additional security when the knob is tightened; the middle jaw locks onto the pins for a more stable construct.

### **About Innovative Medical Products**

Innovative Medical Products Inc.® designs, manufactures, and brings to market patient positioning devices for healthcare and sports related surgeries. With more than 40 years of experience in supporting hospitals, surgery centers, and clinics, IMP provides positioning products that stabilize and protect patients during navigated, computer assisted and robotic surgeries. IMP devices save hospital and surgery center costs by minimizing staff, decreasing liability and improving the overall patient experience. IMP manufactures, patents and trademarks products to Federal and International standards for quality and safety for both the patient and the healthcare institution. For more information, visit: <http://www.impmedical.com>

###

Media contact: Harry McBrien, Maier-True Communication; 860-470-3127; [harry@maier.com](mailto:harry@maier.com)