

DE MAYO

Adapt2Fit[®]

MODULAR KNEE POSITIONER

**Same De Mayo Knee Positioner
with even more versatility.**



90 degree aluminum boot holds patient's heel securely.

Sterile Patient Protective Pads ensures patient comfort and safety.


For use with IMP Universal Distractor.

Boot ball allows for unlimited internal and external rotation.

Longer 32" baseplate provides more flexion and extension travel.



Baseplate separates, allowing for easier cleaning, sterilization, and storage.



KEY FEATURES

- **Enhanced Stability and Visualization**

Creates a more stable operative site with enhanced visualization for all knee replacement surgeries including robotic assisted surgeries

- **More Flexion and Extension**

A 32" baseplate allows for more travel while the leg is supported in the knee positioner

- **Reduced Procedure Fatigue**

The Adapt2Fit Knee Positioner supports the weight of the limb and acts like an extra set of hands, reducing staff fatigue

- **Easier to Clean, Sterilize, and Store**

Knee positioner disassembles into a single container for either blue wrapping or transferring into a rigid container. The smaller footprint also maximizes storage space

ORDERING INFORMATION

Catalog Number	Description
1020*	De Mayo Adapt2Fit Modular Knee Positioner [®]
803-GP-10	Sterile Pressure Protector Pad Set
803-ABD	90° Aluminum Boot for De Mayo Knee Positioner
713-717	De Mayo Push Button™ Pin Locking Clamp
120	Quad Clamp [®]
621	De Mayo Universal Distractor [®] 2.0
907-GP-10	Sterile Pressure Protective Pad Set for De Mayo Universal Distractor [®] 2.0

*Various configurations available. Call IMP for a quote and additional information



IMP[®] products are protected by patent & patent pending rights ~ go to IMPpatents.com
©2023 Innovative Medical Products, Inc.
ALL RIGHTS RESERVED



Positioned to Perform

Innovative
Medical
Products™

87 Spring Lane,
Plainville, CT 06062
PH: (860) 793-0391 F: (866) 459-1805
Toll Free PH: (800) 467-4944
www.IMPmedical.com

SEL-IMP-0020