

MorphBoard®


SURGICAL PEG BOARD

Stable patient positioning in the lateral decubitus position.

PegDok™

PEG BOARD PEG STABILIZER

Pegs are placed into the MorphBoard®, securing the patient in the lateral position.



Center board can be rotated vertically or horizontally to better accommodate patients.



PegDok™ with MorphBoard®



Composite radiolucent pegs are available in sizes ranging from 8" - 14"



MorphBoard® Foam Base Pad has memory foam that improves pressure distribution for the patient.



Tuffease® Gel pad provides an additional layer of protection to prevent pressure ulcers to the patient.



MorphBoard® Leg Extension Kit is used for muscle sparing THA



IMP® products are protected by patent & patent pending rights ~ go to IMPpatents.com
©2025 Innovative Medical Products, Inc.
ALL RIGHTS RESERVED



Positioned to Perform

Innovative
Medical
Products™

87 Spring Lane,
Plainville, CT 06062
PH: (800) 467-4944
www.IMPmedical.com

SEL-IMP-0016 V0

MorphBoard®

SURGICAL PEG BOARD

ORDERING INFORMATION

Catalog Number	Description
105	MorphBoard® Positioning System - 1 ea. Center Board, 2 ea. End Board, 2 ea. Pegs (8", 12", 14"), 2 ea. Clamping Plates and Rail Clamps, 4 PegDoks
105-1	Center Board System - 1 ea. Center Board, 2 ea. Pegs (8", 12", 14"), 2 ea. Clamping Plates and Rail Clamps, 4 PegDoks
105-4	MorphBoard® Leg Extension Kit System - 1 ea. Extension Board, 1 ea. Support Stand, 1 ea. Side Rail, 1 ea. Adjustable Clamping Plate, 1 ea. Rail Clamp
105-2	MorphBoard® End Board 12" x 20"
105-3	MorphBoard® Center Board 30" x 20"
105-5	Leg Extension Board 36" x 10"
105-7	Support Stand w/Leg Extension Board
105-Side Rail	Side Rail w/Leg Extension Board
409	Rail Clamp
805	Adjustable Clamping Plate w/Leg Extension Board
105-8	8" Peg
105-9	9" Peg
105-12	12" Peg
105-14	14" Peg
105-8 CP	8" Composite Peg
105-12 CP	12" Composite Peg
105-14 CP	14" Composite Peg
105-CP	Clamping Plate
105-BP-G	Base Board Pads 18" x 41" x 3/4" (10 per case), Memory Foam
105-6P	6" Peg Foam (20 per case)
105-8P	8" Peg Foam (20 per case)
105-9P	9" Peg Foam (20 per case)
105-12P	12" Peg Foam (20 per case)
105-14P	14" Peg Foam (20 per case)
TE-MBP	TUFFease® Gel Pad for Base Board 50" x 20"
TE-LXP	TUFFease® Gel Pad for Leg Extension 36" x 12"
105-PD-S	PegDok™ Small
105-PD-1	PegDok 1 Hole
105-PD-2	PegDok 2 Hole
105-PD-KIT	PegDok Kit - 2 ea. 105-PD-S, 1 ea. 105-PD-1, 1 ea. 105-PD-2



IMP® products are protected by patent & patent pending rights ~ go to IMPpatents.com
©2025 Innovative Medical Products, Inc.
ALL RIGHTS RESERVED



Positioned to Perform

Innovative
Medical
Products™

87 Spring Lane,
Plainville, CT 06062
PH: (800) 467-4944
www.IMPmedical.com

SEL-IMP-0016 V0