


DE MAYO
Universal Distractor® 2.0
 EXTERNAL JOINT DISTRACTOR

Designed for External Knee Joint Distraction

The De Mayo Universal Distractor® 2.0 applies finite surgeon control pressure to the femur, creating an unobstructed view into the knee joint.



Actual view of joint distraction

Sterile Patient Protective Distractor Pad protects the popliteal area



Hand control for finite adjustment of distraction

Tension release lever to fine-tune spacing or disengage distraction



Multiple positions on the distractor block allow for leg length adjustment

KEY FEATURES

- Compatible with all De Mayo Knee Positioners
- Creates an unobstructed view into the knee joint without the use of lamina spreaders
- Improves access to the posterior aspect of the joint for removal of osteophytes, resected bone, retained cement, and ease of cautery and pain management techniques
- Allows for easy insertion of instruments, trials, and implants while performing Uni, TKA, ACL, and Tibial Plateau surgeries
- Improves minimally invasive surgery techniques for component insertion and prevents component damage
- Create reliable visualization from procedure to procedure for reproducible results

ORDERING INFORMATION

Catalog Number	Description
621	De Mayo Universal Distractor® 2.0 (requires Boot w/Distractor Block)
907-GP-10	Sterile Patient Protective Pad Set for De Mayo Universal Distractor® 2.0



Sterile Patient Protective Pad Set. Latex-free



IMP® products are protected by patent & patent pending rights ~ go to IMPpatents.com
 ©2023 Innovative Medical Products, Inc.
 ALL RIGHTS RESERVED



Positioned to Perform

Innovative Medical Products™

87 Spring Lane,
 Plainville, CT 06062
 PH: (860) 793-0391 F: (866) 459-1805
 Toll Free PH: (800) 467-4944
www.IMPmedical.com