

HUMBLES

LapWrap®

ARM MANAGEMENT SYSTEM

A safe arm management alternative to tucking sheets

Humbles LapWrap is designed to reduce the risks from tucking sheets as described in the AORN perioperative guidelines.



Easy access to all IV lines

Secures arms to O.R. table

Soft & latex free material protects skin

PREVENTS SHEET TUCKING COMPLICATIONS¹

- Hyperextension of the elbow can stretch the median nerve
- If the draw sheet is tucked too tightly, it can create a pressure or nerve injury or a compartment syndrome
- Tucking the patient's arm may interfere with physiological monitoring and could result in an unrecognized infiltrated IV line in the tucked arm
- Tucking the arms has been associated with ulnar neuropathy

ORDERING INFORMATION

Catalog Number	Description
904-G	Humbles LapWrap® Positioning Pad (10/Case)
904-G-EX	Humbles LapWrap® Extensions



Scan to Learn More

KEY FEATURES

- Easily secure patient's arms while allowing accessibility for anesthesiologists
- Adjustable straps allows for easier arm management compared to a continuous draw sheet
- Designed by perioperative professionals to improve patient safety



LapWrap before procedure



LapWrap during procedure

1. AORN Guidelines for Perioperative Practice 2023



IMP® products are protected by patent & patent pending rights ~ go to IMPpatents.com
©2023 Innovative Medical Products, Inc.
ALL RIGHTS RESERVED



Positioned to Perform

Innovative
Medical
Products™

87 Spring Lane,
Plainville, CT 06062
PH: (860) 793-0391 F: (866) 459-1805
Toll Free PH: (800) 467-4944
www.IMPmedical.com

SEL-IMP-0013