

KEEP ALL CRITICAL PARTS TOGETHER

Laser-etched numbers ensure that the correct components (baseplate(s), carriage, boot, and handle) stay together. This extends life, reduces wear, and simplifies assembly of the product.

All critical parts are Serialized with the same unique 3-digit code

1. Carriage
2. Handle
3. Connector Link (only applicable for Adapt2Fit®)
4. Boot
5. Baseplate(s)


Adapt2Fit® (CAT# 1020)



D2® (CAT# 803-114)



V2® (CAT# 114)



De Mayo (CAT# 803)



Catalog Numbers 1020, 803, 803-S, 803-114,
803-114-S, 114, 114-S
Doc Title INFO-IMP-0004
Version 00

IMP® products are protected by patent & patent pending rights ~
go to www.IMPpatents.com
© September 2023 Innovative Medical Products, Inc.
ALL RIGHTS RESERVED