

MORPHBOARD®

LongBoard

SINGLE-PIECE PEGBOARD

Stable patient positioning in the lateral decubitus position.



KEY FEATURES

- **Stable, Single-Piece Design**

The MorphBoard® LongBoard provides a strong, continuous positioning platform without connection points between multiple board sections, helping support a simple set-up before demanding hip procedures.

- **Enhanced Hip Joint Access**

Leg extension facilitates proper lateral patient positioning for minimally invasive or muscle-sparing THA. Improves surgical exposure by enabling controlled posterior leg positioning for anterior hip dislocation.

- **Simplified O.R. Setup**

Integrated clamping areas and a quick-connect design streamline setup on any O.R. table. The single-piece board also eliminates the need to assemble modular sections or separate end-board extensions.



MorphBoard® LongBoard is an alternative that comes in one piece and allows simple setup.



IMP® products are protected by patent
& patent pending rights ~ go to IMPpatents.com
©2026 Innovative Medical Products, Inc.
ALL RIGHTS RESERVED



**Innovative
Medical
Products®**

87 Spring Lane,
Plainville, CT 06062
PH: (800) 467-4944
www.IMPmedical.com

SEL-IMP-0023 R00

MORPHBOARD®

LongBoard

SINGLE-PIECE PEGBOARD

ORDERING INFORMATION

Catalog Number	Description
105-LB	MorphBoard® LongBoard
105-LB KIT	MorphBoard® LongBoard System - 1 LongBoard, 6 Composite Pegs, 2 Drop Bars, 2 Rail Clamps, 1 PD-KIT
105-LB-EXT	MorphBoard® LongBoard Leg Extension
105-LB-EXT KIT	MorphBoard® LongBoard System - 1 LongBoard Leg Extension, 1 Rail Clamp, 1 Adjustable Clamping Plate, 1 LongBoard Leg Extension Hinge Pin
105-LB-EXT-HP	MorphBoard® LongBoard Leg Extension Hinge Pin
409	Rail Clamp
805	Adjustable Clamping Plate
105-8	8" Peg
105-9	9" Peg
105-12	12" Peg
105-14	14" Peg
105-8 CP	8" Composite Peg
105-9 CP	9" Composite Peg
105-12 CP	12" Composite Peg
105-14 CP	14" Composite Peg
105-CP	Clamping Plate
105-BP-G	Base Board Pads 18" x 41" x 3/4" (10 per case), Memory Foam
105-6P	6" Patient Protective Peg Pad (20 per case)
105-8P	8" Patient Protective Peg Pad (20 per case)
105-9P	9" Patient Protective Peg Pad (20 per case)
105-12P	12" Patient Protective Peg Pad (20 per case)
105-14P	14" Patient Protective Peg Pad (20 per case)
TE-MBP	TUFFease® Gel Pad for Base Board 50" x 20"
TE-LXP	TUFFease® Gel Pad for Leg Extension 36" x 12"
105-PD-S	PegDok™ Small
105-PD-1	PegDok 1 Hole
105-PD-2	PegDok 2 Hole
105-PD KIT	PegDok Kit - 2 ea. 105-PD-S, 1 ea. 105-PD-1, 1 ea. 105-PD-2



IMP® products are protected by patent
& patent pending rights ~ go to IMPpatents.com
©2026 Innovative Medical Products, Inc.
ALL RIGHTS RESERVED



Positioned to Perform

Innovative
Medical
Products®

87 Spring Lane,
Plainville, CT 06062
PH: (800) 467-4944
www.IMPmedical.com

SEL-IMP-0023 R00