

MorphBoard®

SURGICAL PEG BOARD

Stable patient positioning in the lateral decubitus position.

PegDok™

PEG BOARD PEG STABILIZER

KEY FEATURES

• Optimal Patient Positioning

The MorphBoard®'s modular components allow simple conversion for bariatric patient support on any O.R. table, while multiple peg lengths provide greater versatility

• Enhanced Peg Stability

New PegDok™ reduces peg dislodging with enhanced peg stability


• The Widest Range of Patients

Along with its two modular end boards, the MorphBoard® surgical peg board features a lightweight center board that can be easily adjusted to accommodate a bariatric patients

• Easy Setup and Storage

Lightweight, modular components make set up and storage simple

Pegs are placed into the MorphBoard®, securing the patient in the lateral position.



Center board can be rotated vertically or horizontally to better accommodate patients.



PegDok™ with MorphBoard®



Composite radiolucent pegs are available in sizes ranging from 8"- 14"

MorphBoard® Foam Base Pad has memory foam that improves pressure distribution for the patient.



Tuffease® Gel pad provides an additional layer of protection to prevent pressure ulcers to the patient.



MorphBoard® Leg Extension Kit is used for muscle sparing THA



IMP® products are protected by patent & patent pending rights ~ go to IMPpatents.com
©2026 Innovative Medical Products, Inc.
ALL RIGHTS RESERVED

IMP
Positioned to Perform

Innovative Medical Products®

87 Spring Lane,
Plainville, CT 06062
PH: (800) 467-4944
www.IMPmedical.com

SEL-IMP-0016 R01

MorphBoard®

SURGICAL PEG BOARD

ORDERING INFORMATION

Catalog Number	Description
105	MorphBoard® Positioning System - 1 ea. Center Board, 2 ea. End Board, 2 ea. Pegs (8", 12", 14"), 2 ea. Clamping Plates and Rail Clamps, 4 PegDoks
105-1	Center Board System - 1 ea. Center Board, 2 ea. Pegs (8", 12", 14"), 2 ea. Clamping Plates and Rail Clamps, 4 PegDoks
105-4	MorphBoard® Leg Extension Kit System - 1 ea. Extension Board, 1 ea. Support Stand, 1 ea. Side Rail, 1 ea. Adjustable Clamping Plate, 1 ea. Rail Clamp
105-2	MorphBoard® End Board 12" x 20"
105-3	MorphBoard® Center Board 30" x 20"
105-5	Leg Extension Board 36" x 10"
105-7	Support Stand w/Leg Extension Board
105-Side Rail	Side Rail w/Leg Extension Board
409	Rail Clamp
805	Adjustable Clamping Plate w/Leg Extension Board
105-8 CP	8" Composite Peg
105-12 CP	12" Composite Peg
105-14 CP	14" Composite Peg
105-CP	Clamping Plate
105-BP-G	Base Board Pads 18" x 41" x 3/4" (10 per case), Memory Foam
105-6P	6" Peg Foam (20 per case)
105-8P	8" Peg Foam (20 per case)
105-9P	9" Peg Foam (20 per case)
105-12P	12" Peg Foam (20 per case)
105-14P	14" Peg Foam (20 per case)
TE-MBP	TUFFease® Gel Pad for Base Board 50" x 20"
TE-LXP	TUFFease® Gel Pad for Leg Extension 36" x 12"
105-PD-S	PegDok™ Small
105-PD-1	PegDok 1 Hole
105-PD-2	PegDok 2 Hole
105-PD-KIT	PegDok Kit - 2 ea. 105-PD-S, 1 ea. 105-PD-1, 1 ea. 105-PD-2



IMP® products are protected by patent & patent pending rights ~ go to IMPpatents.com
©2026 Innovative Medical Products, Inc.
ALL RIGHTS RESERVED



Positioned to Perform

Innovative
Medical
Products®

87 Spring Lane,
Plainville, CT 06062
PH: (800) 467-4944
www.IMPmedical.com

SEL-IMP-0016 R01